

## Day-1

### DAY 1 (MON DEC 16): GET IT WORKING

### BEFORE HACKATHON STARTS (Do this Sunday night!)

#### Setup Checklist:

```
# 1. Create project folder
mkdir MicroscoPilot
cd MicroscoPilot

# 2. Create Python virtual environment
python -m venv venv
source venv/bin/activate # On Windows: venv\Scripts\activate

# 3. Install basics
pip install jupyter numpy matplotlib anthropic pillow
```

#### Get Anthropic API Key:

1. Go to: <https://console.anthropic.com/>
2. Sign up (you get \$5 free credits!)
3. Create API key
4. Save it in `.env` file:

```
ANTHROPIC_API_KEY=your_key_here
```

#### Clone DTMicroscope:

```
git clone https://github.com/pycroscopy/DTMicroscope.git
cd DTMicroscope
pip install -e .
cd ..
```

#### Test it works:

```
cd DTMicroscope/notebooks
jupyter notebook
```

```
# Open and run: 01-AFM-Demo.ipynb
# If you see a microscope image, you're good!
```

---

## **Hour 1-3 (9:00 AM - 12:00 PM): Opening & Foundation**

### **9:00-10:30 AM: Watch Opening Ceremony**

- Take notes on challenges they mention
- Screenshot any specific tasks
- Note judging criteria
- Join Slack immediately

### **10:30 AM - 12:00 PM: Build Project Skeleton**

#### **PROMPT #1 TO CLAUDE/CURSOR:**

I'm building MicroscoPilot for the Microscopy Hackathon. I'm a beginner in ML.

PROJECT: Autonomous microscopy agent using Claude Vision API + DTMicroscope simulator

#### **REQUIREMENTS:**

- Use DTMicroscope's AFM digital twin (already installed)
- Agent makes decisions using Claude Sonnet 4 vision API
- Simple autonomous exploration: move around, capture images, find interesting features
- Track discoveries and show visualizations

#### **CREATE:**

Project structure:

/MicroscoPilot

  /src

    \_\_init\_\_.py

    agent.py           # Main agent logic

    vision.py          # Claude vision API calls

    microscope.py      # DTMicroscope wrapper

    memory.py          # Store what agent has seen

    visualizer.py      # Display results

  /tests

```
test_basic.py
/outputs      # Save results here
/data         # Sample images
main.py       # Run the agent
requirements.txt
README.md
.env.example
.gitignore
```

Use SIMPLE Python (I'm learning!):

- Clear variable names
- Lots of comments explaining what each part does
- Error handling with helpful messages
- Logging to see what's happening

Start with main.py that imports everything and has a basic run loop.

**WHAT YOU'LL GET:** Full project structure. Test that everything imports correctly.

---

## HOURL 4-6 (12:00 PM - 2:00 PM): VISION MODULE

### PROMPT #2:

Build src/vision.py – the module that uses Claude Vision API to analyze microscope images.

#### REQUIREMENTS:

1. Function: analyze\_image(image\_path) -> dict
  - Loads image from file
  - Converts to base64 (required for API)
  - Sends to Claude Sonnet 4 with specialized prompt
  - Returns structured analysis
2. Prompt engineering (CRITICAL – this is where intelligence comes from):
  - Tell Claude it's analyzing Atomic Force Microscopy (AFM) images and should:
    - Identify features: grain boundaries, defects, particles, rough regions, smooth regions
    - Rate image quality (0-10)

- Assess interestingness (is there something worth investigating?)
- Suggest next action (zoom in, move left/right/up/down, adjust settings)
- Provide confidence score
- Explain reasoning

### 3. Response parsing:

- Extract Claude's response into Python dict
- Handle errors gracefully (API might be slow or fail)
- Log the analysis to a file

### 4. Include retry logic (API calls can fail):

- Try up to 3 times
- Wait 2 seconds between retries
- Give helpful error message if all fail

### TEACHING NOTES:

- Explain what base64 encoding is (turning image into text)
- Show example of the API call format
- Add comments explaining each step
- Include a test function at the bottom that analyzes a sample image

IMPORTANT: Make it work with FREE API tier (don't spam requests)

**WHAT YOU'LL GET:** Working vision module. Test with:

```
python src/vision.py --test data/sample_afm_image.png
```

---

## **HOUR 7-8 (2:00 PM - 4:00 PM): MICROSCOPE CONTROLLER**

### PROMPT #3:

Build src/microscope.py – wrapper around DTMicroscope for easy control.

### REQUIREMENTS:

#### 1. Class: MicroscopeController

- Initialize DTMicroscope AFM simulator
- Provide simple methods like:
  - \* move\_to(x, y) – move stage to coordinates

- \* `capture_image()` – get current image
- \* `get_current_position()` – where are we?
- \* `zoom_in()` / `zoom_out()` – adjust magnification
- \* `get_state()` – full microscope status

2. Safety checks:

- Don't allow moving outside valid boundaries (0 to `sample_size`)
- Validate all inputs
- Return clear error messages

3. State tracking:

- Remember where we've been (position history)
- Track number of images captured
- Store settings used for each image

4. Simulation mode (for testing without `DTMicroscope`):

- If `DTMicroscope` not available, create fake microscope
- Generate random but realistic images
- Still track positions correctly

TEACHING NOTES:

- Explain what a "wrapper class" is (simpler interface to complex code)
- Show how to check if a library is installed
- Document the coordinate system (x,y axes)
- Add diagram in comments showing how positions map to sample

MAKE IT BEGINNER-FRIENDLY:

- Use descriptive method names
- Include docstrings with examples
- Add assertions to catch mistakes
- Print warnings for weird inputs

**WHAT YOU'LL GET:** Microscope control module. Test with:

```
python src/microscope.py --demo
```

---

## HOURL 9-10 (4:00 PM - 6:00 PM): BASIC AGENT BRAIN

### PROMPT #4:

Build src/agent.py – the autonomous decision-making agent.

This is THE CORE of your project. The agent needs to:

1. Explore the sample intelligently
2. Recognize interesting features
3. Learn from what it sees
4. Make decisions about where to go next

REQUIREMENTS:

Class: AutonomousAgent

```
__init__(microscope, vision_analyzer):  
    - Takes microscope controller and vision module  
    - Initialize empty memory  
    - Set starting strategy (e.g., random exploration)
```

```
explore(num_steps=100):  
    - Main loop that runs autonomously  
    - For each step:  
        1. Capture current image  
        2. Analyze with vision AI  
        3. Decide next action based on analysis  
        4. Execute action  
        5. Store results in memory  
        6. Print progress update
```

```
decide_next_action(vision_analysis, current_state):  
    - Decision logic:  
        * If interesting (score > 7): zoom in and capture more  
        * If seen before: skip and move elsewhere  
        * If quality is poor: adjust settings  
        * If nothing interesting for N steps: jump to random new area  
        * If near boundary: turn around  
    - Return action: {'type': 'move', 'x': 10, 'y': 20}
```

TEACHING MOMENT – Explain the agent loop:

"The agent is just a while loop! It keeps running until you tell it to stop. Each loop: sense the world (vision) → think (decide\_next\_action) → act (move microscope)"

KEEP IT SIMPLE FOR NOW:

- No fancy ML yet, just if/else logic
- Random exploration is fine
- Focus on getting the loop working
- We'll add intelligence tomorrow

#### LOGGING:

- Print emoji-based updates so demos look cool:
  - "🔍 Analyzing position (45, 120)..."
  - "✨ Interesting feature detected! Moving closer..."
  - "📍 Exploring new region..."
  - "😴 Nothing interesting here, moving on..."

**WHAT YOU'LL GET:** Basic working agent! Test with:

```
python main.py --steps 20
```

---

## HOOR 11 (6:00 PM - 7:00 PM): INTEGRATION TEST

### Your Tasks (No Claude needed):

#### 1. Create main.py:

```
# This is where everything comes together
from src.microscope import MicroscopeController
from src.vision import VisionAnalyzer
from src.agent import AutonomousAgent

# Initialize components
microscope = MicroscopeController()
vision = VisionAnalyzer(api_key="your_key")
agent = AutonomousAgent(microscope, vision)

# Run exploration
print("🚀 Starting autonomous exploration...")
results = agent.explore(num_steps=50)
print(f"✅ Complete! Found {len(results['discoveries'])} interesting features")
```

#### 2. Test everything:

```
python main.py
```

**EXPECTED:** Agent moves around, captures images, makes decisions. You see progress logs.  
No crashes!

3. **Debug if needed** (save Claude messages for actual bugs):

- Read error messages carefully
- Check that all imports work
- Verify API key is loaded
- Make sure DTMicroscope is accessible

---

## Hour 12 (7:00 PM - 8:00 PM): DINNER + TEAM CHECK-IN

### Your Tasks:

- EAT REAL FOOD (seriously, brain fuel!)
- Review what works and what doesn't
- Push code to GitHub
- Post in Slack: "Built autonomous microscopy agent. Looking for teammates interested in [ML/microscopy/visualization]. Have working prototype!"
- Plan tomorrow's features

### GitHub Setup:

```
git init
git add .
git commit -m "Day 1: Basic agent working"
git remote add origin https://github.com/yourusername/MicroscoPilot.git
git push -u origin main
```

---

## Hour 13-14 (8:00 PM - 10:00 PM): BASIC VISUALIZATION

### PROMPT #5:



Build `src/visualizer.py` – real-time visualization of agent exploration.

## REQUIREMENTS:

Create a dashboard using `matplotlib` that shows:

### 1. Main window with 4 panels:

TOP LEFT: Current microscope view

- Show the latest image agent captured
- Mark current position

TOP RIGHT: Exploration map

- 2D heatmap of where agent has been
- Color: blue = never visited, yellow = visited once, red = visited many times
- Show current position as red dot
- Mark discovered features with stars 🌟

BOTTOM LEFT: Discoveries timeline

- Line graph: time vs. number of features found
- Show when interesting things were discovered

BOTTOM RIGHT: Statistics

- Text display showing:
  - \* Total steps taken
  - \* Areas explored
  - \* Features found
  - \* Current strategy
  - \* Time elapsed

### 2. Real-time updates:

- Update every time agent takes a step
- Don't freeze (use `matplotlib.animation` or just `plt.pause()`)
- Save final state as PNG

### 3. Export capability:

- Save exploration video as MP4 (optional, might need `ffmpeg`)
- Save heatmap as high-res image
- Export data as JSON

## TEACHING NOTE:

"Visualization is CRITICAL for demos! Judges need to SEE your agent working. This is the difference between 'cool code' and 'OH WOW' reaction."

KEEP IT FUNCTIONAL, NOT FANCY (for now):

- Basic matplotlib is fine
- Don't worry about making it beautiful yet
- Focus on showing the agent's behavior clearly
- Tomorrow we'll make it pretty

**WHAT YOU'LL GET:** Live visualization. Test with:

```
python src/visualizer.py --demo # Shows fake data
python main.py --visualize      # Shows real agent
```

---

## END OF DAY 1 CHECKLIST

- ✓ ~~Project structure created~~
- ✓ ~~Vision module working (can analyze images with Claude)~~
- ✓ ~~Microscope controller working (can move and capture)~~
- ✓ ~~Basic agent working (explores autonomously)~~
- ✓ ~~Visualization shows agent behavior~~
- ✓ ~~Code on GitHub~~
- ✓ ~~You understand what each part does~~
- ✓ ~~You have a demo you can show~~

**EXPECTED STATE:** You have a working autonomous agent that:

- Moves around the microscope sample
- Captures images
- Uses Claude Vision to analyze them
- Makes simple decisions
- Shows its exploration visually

**It's basic, but it WORKS!** That's more than many teams will have.

---